

Hoof, skin & wounds

Sanocyn forte Hydrogel

Sanocyn forte Hydrogel forms a protective film on injured skin and keeps it moist. The gel is preserved by hypochlorous acid (HOCl). Its adhesive consistency means it stays in place for a long time, even in hard-to-reach areas. It complements Sanocyn forte Solution as a system for skin, hoof and claw.

Pack sizes

75 ml | **350 ml**

Composition	**Water, sodium chloride, active chlorine released from hypochlorous acid (HOCl), layered silicate (LiMgNa silicate). Produced in a proprietary electrochemical process**
Active chlorine from HOCl	**0.5 g/l (< 0.1 % w/v)**
Preservation	**by hypochlorous acid (HOCl)**
pH	**approx. 6.5 (near-neutral)**

Skin compatibility	**Sanocyn forte Solution and Hydrogel were classified as non-irritant to skin in the in vitro skin irritation test according to ISO 10993-10 / OECD 439 (reconstructed human epidermis; TentaMedix GmbH, DAkkS-accredited, 2019).**
Shelf life	**unopened 24 months, 3 months after opening (according to stability data)**
	Free from alcohol, antibiotics, colourings and fragrances
Species	**all mammals of any age (according to the product information)**
	Manufactured and filled in Austria and Germany
EAN	**75 ml 9120092250020 · 350 ml 9120092250051**

AREAS OF USE

Areas of use

01 Protective film on injured skin

02 Keeps the area moist over a longer period

03 Care of hard-to-reach areas

04 Care of hoof and frog

05 Keeping wounds moist (moist wound environment), can be combined with standard wound dressings

DIRECTIONS FOR USE

Directions for use

Recommended combination

> Recommended as a system: rinse with Sanocyn forte Solution first, then apply the Hydrogel after one to two minutes.

01 Expose the area and clean it with Sanocyn forte Solution. The gel adheres better to dry skin.

02 Shake the bottle before use.

03 Apply the Hydrogel in a targeted way and generously. It is absorbed slowly and does not need to be rinsed off.

04 Cover with a bandage if needed. If contamination is not expected, the area can be left uncovered.

05 Repeat several times a day.

01 Clean the hoof or claw carefully with a hoof pick.

02 Shake the bottle before use.

03 Spray hoof and frog generously with the Hydrogel.

04 Afterwards, keep the animal on clean, dry bedding where possible.

05 Repeat several times a day.

Notes

- Follow the veterinarian's instructions where applicable.
- Read the labelling and product information before use.
- For external (topical) use only. In case of eye contact, rinse with plenty of water.
- Store at room temperature and protect from direct sunlight.
- Do not decant into other containers.
- Do not use in animals with known hypersensitivity to the ingredients.
- Shake the bottle before each use.
- Do not use together with other wound-cleansing products or with solutions containing silver, iodine or polyhexanide. Active chlorine has an oxidising effect.