

Professional Use

Sanocyn Pro

Sanocyn Pro is the 500 ml container for professional use in the practice. The wound rinsing solution is preserved by hypochlorous acid (HOCl). It is intended for cleansing and rinsing wounds in the practice, not for rinsing body cavities.

Pack sizes

500 ml, Professional Use

Composition	**Water, sodium chloride, active chlorine released from hypochlorous acid (HOCl)**
Active chlorine from HOCl	**0.5 g/l**
pH	**near-neutral**
Shelf life	**unopened 24 months, 3 months after opening (according to stability data)**
	Free from antibiotics and alcohol
	Manufactured and filled in Austria and Germany

EAN	9120092250853

AREAS OF USE

Areas of use

01 Wound rinsing solution in a 500 ml practice container for cleansing and rinsing wounds, not for rinsing body cavities

02 Rinsing of hoof and claw

DIRECTIONS FOR USE

Directions for use

Recommended combination

For follow-up care at home, the practice can supply Sanocyn forte Solution and Hydrogel. Sanocyn Pro can be combined with Sanocyn forte Hydrogel.

01 Expose and clean the area.

02 Rinse generously with Sanocyn Pro.

03 Cover with suitable dressing material if needed.

04 To make wound dressings easier to remove, soak them with Sanocyn Pro.

Notes

- Professional Use.
- Sold through medical professionals only.

- Follow the veterinarian's instructions where applicable.
- Store at room temperature and protect from direct sunlight.
- For external use only.
- Dispose of empty containers, fully emptied, with household waste.